



OWL *On Wall Logic* Touch™

Columbia Dentoform® Teaching Solutions
DentalEZ® Equipment
NevinLabs™ Workstations
RAMVAC® Utility
StarDental® Instruments



The OWL Touch™ incorporates a new full color, touch screen display to provide easy control and monitoring of your dental utility room equipment. Simple, easy to understand screen layouts make navigation of this powerful device a cinch.

**Demo the OWL Touch experience
online at DentalEZ.com/OWL.**

RAMVAC® Utility Solutions

Four different categories of devices can be controlled through a single OWL Touch panel:

Vacuum Units

RAMVAC® brand vacuum units with the new S2 control can be easily plugged into the OWL Touch and will display maintenance notifications, fault messages and other system data. However, even non-RAMVAC vacuum units can be controlled via an OWL Interface Control attached to the OWL Touch.

Air Compressors

RAMVAC Osprey™ air compressors featuring the new C2 control (Smart control) are easily connected and will provide similar information to the vacuum units. Osprey compressors with a Basic control or non-RAMVAC brand compressors can also be connected via the OWL Interface Control.

Water Valve

Any type of water valve utilizing a solenoid valve can be easily turned off and on utilizing an OWL Interface Control.

Auxiliary Devices

Any type of system that can be turned off or on with a switch can be controlled by an OWL Interface Control (up to 10 amps). Systems such as lighting, stereo or intercom can be attached to the OWL Touch and turned on and off either manually or with the schedule.



RAMVAC products are manufactured in an ISO 13485:2003 Certified facility.

☎ 866.DTE.INFO

🌐 www.DentalEZ.com

📺 Follow Us! 🐦



OWL *On Wall Logic* Touch™

- Home screen provides the ability to choose Manual or Schedule run mode for the equipment, and includes a Master Shutdown button.
- Manual mode allows the user to turn on or off each individual type of device, similar to a standard remote switch panel.
- Schedule mode automatically turns on or off all devices by a saved user defined weekly schedule.
- Master Shutdown turns off the entire utility room when activated. Simply push this button and your utility room stays off until either the Manual or Schedule mode is selected.
- OWL Touch models are available with either a four, eight or sixteen port OWL Hub device.
- Entire configuration from device, to hub, to OWL is easily connected with category 6 patch cable.



OWL configurations:

OWL Touch models are available with 4, 8, or 16 port OWL Hubs.

